

An underwater scene featuring a vibrant coral reef in the foreground and numerous small, colorful fish swimming in the clear blue water above. The image is overlaid with a semi-transparent blue gradient.

Microplastics Analysis


— Analytical methods with FTIR and Raman —

Microplastics Analysis

Microplastics, which are generally less than 5 mm in size, have raised concerns about accidental ingestion by marine organisms, contamination of commercial drinking water, and health hazards to those who consume them. As a result, there is a growing movement to investigate the effects of the microplastics. JASCO FTIR and Raman spectrometers can easily perform various types of microplastics analysis.

Use of different instruments depending on size


JASCO offers instruments capable of analyzing microplastics in a wide range of sizes from a few millimeters to sub-micron orders.




FTIR spectrometer

FTIR


FTIR spectrometer can be used for qualitative analysis of microplastics over several hundred microns. The measurement results of polyethylene (PE) and polyethylene terephthalate (PET) are shown here.




FTIR microscope can be used to analyze microplastics as small as 5 μm to 10 μm in length. Qualitative analysis is of course possible, imaging measurements can be used to visualize the distribution of components, and particle size analysis is also possible. These are example imaging analysis data about several kinds of filtrated samples using of polytetrafluoroethylene (PTFE) filter as polyethylene (PE), polypropylene (PP), polystyrene (PS), polyethylene terephthalate (PET) and polyvinyl chloride (PVC) in water.



Averaged spectra of each component



Chemical Image (peak height)





Stacked histograms by component

Those samples and data obtained by FTIR microscope were provided by TOSOH Analysis and Research Center Co., Ltd.

Application Data

100-MT-0262

Raman microscope can be used for microplastics sample analysis as sub-micron order. In addition to qualitative analysis, mapping measurement enables visualization of component distribution as well as FTIR microscope, and particle size analysis is also available. Here are the results of our measurements of Airborne Microplastics (AMPs).



Measurement Results of Airborne Microplastics (AMPs)

Those samples were provided by Prof. Okochi of Waseda University.

Application Data

100-AN-0031

Macro analysis of Microplastics

Microplastics from several mm to several hundred μm can be easily measured using FTIR ATR method. ATR is a contact sampling method in FTIR reflection measurement.

Fourier Transform Infrared Spectrometer : FTIR

The FT/IR-4X is a powerful Mid-IR FTIR spectrometer, with many features that you find in a research grade instrument, which has high resolution, High S/N and high sensitivity analysis performance.



FT/IR-4X with ATR PRO 4X VIEW

ATR PRO 4X VIEW : ATR PRO ONE X VIEW

ATR PRO 4X VIEW / ATR PRO ONE X VIEW are ATR accessories for FTIR that allow you to measure sample while observing sample surface through crystal. This observation function helps to confirm proper sample contact condition even for several hundred μm order samples and enables to measure it. Images and spectra of samples closely adhering to the ATR crystal can be saved in a single file.



ATR PRO 4X VIEW and
ATR PRO ONE X VIEW



Close contact with
ATR crystal sample (PP) image

FTIR

Application Data

030-AT-0275

USB Microscope kit

This USB Microscope kit is a complete package for microscope observation for FT/IR-4X. This kit can be added to ATR accessory and observation view on the top of crystal can be obtained before measurement. Furthermore, such observation view data and spectra data can be saved in the same file. Even if sample is quite small, this kit supports to handle sample setting with each because sample can be set on ATR crystal while checking the observation view from USB microscope.

* Any accessories from ATR PRO 4X, ATR PRO 4X VIEW, ATR PRO ONE X and ATR PRO ONE X VIEW are necessary.



FT/IR-4X with USB Microscope kit

FTIR

SWITCH with EASE
between

USB
microscope

IMAGE
ACQUISITION

&

ATR

MEASUREMENT

USB microscope
IMAGE
ACQUISITION



ATR
MEASUREMENT




USB microscope image

Measurement Spectrum

ATR close-up image*

SAVE
in 1 FILE



Sample: Microplastics collected from lake water


* Sample observation image that contacted with ATR crystal can be obtained with use of ATR PRO 4X VIEW or ATR PRO ONE X VIEW.

Those samples were provided by e-kagaku Biwako Environment Project 2021 (SEP2021).

Application Data

030-AT-0275

It has two functions "Classification" and "Spectra Search", which support analyzing spectra easily and precisely. The "Classification" function recognizes 35 types of compounds based on a machine learning from various materials and high molecular spectra. It also shows key bands of compound so that an operator confirm a validity of results. The "Spectra Search" function performs searching from approximately 600 spectra library.



1 Measured spectrum

Easy to compare with 2

2 Reference spectrum

A reference spectrum of the selected category

3 Result list

Score and sample name are displayed

4 Advice function

Supports for understanding the result

Some microplastics are degraded in a certain environment and it might appear on the spectrum changes. It may show incorrect search results on a spectrum library search. The "Classification" function still works in such conditions because it focuses on key bands on the spectrum and is not affected by the degradation. The Advanced Spectra Search is a useful software for analyzing microplastics in various types and conditions.

Application Data

280-AN-0274

IR plastics library

This is a dedicated library for microplastics including 39 types, 153 IR spectra. It improves the search result of a spectra library. It also includes cellulose and rubbers in the library, which are typically found in environmental conditions to identify whether it is plastic or not.

Library list

Low-density polyethylene, LDPE	High-density polyethylene, HDPE	Polypropylene, PP	Polyvinyl chloride, PVC
Polystyrene, PS	Foamed polystyrene	Polyamide, PA	Polycarbonate, PC
Poly(ethylene terephthalate), PET	Polyester	Polyurethane, PU	Poly(methyl methacrylate), PMMA
Polyvinyl alcohol, PVA	Polyimide	Poly(butylene succinate)-poly(lactic acid) copolymer, PBS-PLA copolymer	Artificial turf
Polytetrafluoroethylene, PTFE	Acrylonitrile butadiene styrene, ABS	Epoxy resin, EP	Phenol formaldehyde, PF
Polyethersulfone, PES	Polyetherimide, PEI	Polyphenylene sulfide, PPS	Polybutylene terephthalate, PBT
Acrylic fiber	Cellulose acetate	Aramid	Cupro, Cuprammonium rayon
Polynosic, Modal	Viscose rayon	Silk	Wool
Cotton	Cellulose	Styrene-butadiene rubber, SBR	Ethylene-propylene rubber
Isobutylene-isoprene rubber, IIR	cis-1,4-Polyisoprene rubber	trans-1,4-Polyisoprene rubber	

Application Data

280-AN-0274

Microplastics obtained with filtration

Microplastics smaller than hundreds of μm can be measured with the micro FTIR and micro Raman while they are placed on the membrane filter. JASCO micro FTIR and Raman measure one/multi points, imaging with the exactly same point, area between both instruments.

Micro FTIR spectrometer

Micro FTIR observes and measures small samples with microscope function. It also performs imaging measurements, creates chemical images and particle analysis. Micro FTIR obtains information of materials such as polyethylene and polypropylene.



IRT-7200 with FT/IR-4X

Micro Raman spectrometer

Micro Raman has higher spatial resolution compared to the micro FTIR, so it can measure microplastics in the μm range which are difficult to measure with the micro FTIR. The micro Raman also detects additives, pigments in addition to base materials.



NRS-4500

VFS-25-1 Filtration set

FTIR

IRT

NRS

Combined with the filtration set and 25 mm dia. membrane filter, microplastics are filtered on it. Various types of materials and pore size filters are available (sold separately).

* Aspirator or dry scroll pump is required separately.

Lineup of membrane filter

No.	Material	Pore size
1	P T F E	10 μm
2	P T F E	5 μm
3	P T F E	1 μm
4	Alumina	0.2 μm



Membrane filter
(sold separately)

Filtration set

Application Data

110-MT-0271


SH02-FH Filter holder

IRT

NRS

This is a dedicated holder for the 25 mm dia. membrane filter. It can be set on the sharing holder and used on the micro FTIR and Raman.

* SH-02 sharing holder plate or IQ frame is required.



Set the membrane filter.



Fix with the plate cover.



Set on the sharing holder.

Application Data

030-MT-0276

030-AN-0040

Conditate measurement positions are automatically recognized in real time. These dedicated points can be registered as measurement positions at once. It is possible to target measurement positions in detail by setting parameters related to such characteristics. The acquired spectrum is matched with JASCO's original database, and the names of the hit substances are displayed on the observation image.

* Only the NRS series can narrow down the number of measurements point.

NARROW DOWN THE MEASUREMENT POINTS

TARGETING by SIZE

TARGETING by COLOR

AUTOMATIC RECOGNITION

Real-time component identification


Application Data: 260-AN-0018

The iQ Frame easily measures the same minute region and the same point IR spectroscopy and Raman spectroscopy. Both methods are known as a complementary method to analyze information of molecules. For this reason, the iQ Frame has a multilateral analysis result can be obtained by combining both measurements.



IR spectroscopy


Raman spectroscopy



Example of composite analysis of IR and Raman by the iQ Frame.

Application Data: 100-MT-0270, 100-AN-0034


The JASCO Particle Analysis program enables the acquisition of shape and color information about sample from color maps and observation images. Histograms, frequency distributions, and correlation distributions can be created based on the obtained information. Additionally, it is also possible to create the stacked histograms for each component.



Measurement result

No.	Image	Centroid-X	Centroid-Y	Size
1		9.8	8.2	10.482
2		16.2	15.2	9.480
3		102.6	17.3	9.364
4		26.7	30.7	9.667
5		39.4	40.4	9.648
6		42.1	50.0	9.641

Character of each particles



Stacked histogram

Application Data: 100-MT-0270, 100-AN-0034

LIST OF ACCESSORIES FOR EACH MODEL

FT/IR-X series (FT/IR-4X, FT/IR-6X, FT/IR-8X)

ATR PRO 4X	ATR PRO 4X Single reflection ATR accessory with iQX Accessory, for FT/IR-4X only
ATR PRO 4X VIEW	ATR PRO 4X VIEW Single reflection ATR accessory with camera with iQX Accessory, for FT/IR-4X only
ATR PRO ONE X	ATR PRO ONE X Single reflection ATR accessory with iQX Accessory, for 6X/8X
ATR PRO ONE X VIEW	ATR PRO ONE X VIEW Single reflection ATR accessory with camera with iQX Accessory, for 6X/8X
	USB Microscope set (Custom-made) Components <ul style="list-style-type: none"> USB Microscope for observation USB Microscope-ATR Mounting unit for ATR PRO 4X, ATR PRO 4X VIEW, ATR PRO ONE Mask for microscope observation Micro Imaging Analysis Program (including Scale Correction function) * ATR PRO 4X, ATR PRO 4X VIEW, ATR PRO ONE or ATR PRO ONE VIEW is required.
ADSS-4X	ADSS-4X Advanced Spectra Search program for Spectra Manager 2.5
LIB-PLA-ADSS	LIB-PLA-ADSS IR plastics library * ADSS-4X is required

IRT series and NRS series Microscopes

JPA	JASCO Particle Analysis
SH02-FH	SH02-FH Filter holder * SH02 Sharing holder plate (SH02, SH02-IR or SH02-RM) is required.
SH02	SH02 Sharing holder plate (plate only)
VFS-25-1	VFS-25-1 Filtration set for micro plastic preparation, 25 mm dia. Components <ul style="list-style-type: none"> Suction filter bottle, Funnel, Base, Clamp, Silicone stopper * This set does not include an aspirator or dry vacuum pump.
MF25-PTFE10	MF25-PTFE10, PTFE Membrane Filters for VFS-25
MF25-PTFE5	MF25-PTFE5, PTFE Membrane Filters for VFS-25
MF25-PTFE1	MF25-PTFE1, PTFE Membrane Filters for VFS-25
MF25-AL02	MF25-AL02, Aluminum Oxide Membrane Filters for VFS-25

IRT series (IRT-7200, IRT-7100, IRT-5200)

ASN-IR-X	ASN-IR-X IR Advanced Search NAV for Spectra Manager 2.5
ADSS-4X	ADSS-4X Advanced Spectra Search program for Spectra Manager 2.5
LIB-PLA-ADSS	LIB-PLA-ADSS IR plastics library * ADSS-4X is required.
SH02-IR	Sharing holder & iQ Frame for IR

NRS series (NRS-4500, NRS-5500, NRS-7500)

ASN-457	Advanced Search NAV
SH02-RM	Sharing holder & iQ Frame for Raman



JASCO CORPORATION

2967-5, Ishikawa-machi, Hachioji-shi, Tokyo 192-8537 Japan
 Tel: +81-42-649-5177 Fax: +81-42-646-4515 Web: www.jasco.co.jp
 Japan

JASCO INTERNATIONAL CO., LTD.

11-10, Myojin-cho 1-chome, Hachioji-shi, Tokyo 192-0046, Japan
 Tel: +81-42-649-3247 Fax: +81-42-649-3518 Web: www.jascoint.co.jp/english/
 Australia, Hong Kong, India, Indonesia, Korea, Malaysia, New Zealand, Pakistan,
 Philippines, Russia and CIS countries, Singapore, Taiwan, Thailand, Vietnam

JASCO INCORPORATED

28600 Mary's Court, Easton, Maryland 21601, U.S.A.
 Tel: +1-410-822-1220 Fax: +1-410-822-7526 Web: www.jascoinc.com
 Argentina, Bolivia, Brazil, Canada, Chile, Colombia, Costa Rica, Guatemala, Mexico,
 Paraguay, Peru, Puerto Rico, United States of America, Uruguay, Venezuela

JASCO EUROPE S.R.L.

Via Luigi Cadorna 1, 23894 Cremella (LC), Italy
 Tel: +39-039-9215811 Fax: +39-039-9215835 Web: www.jascoeurope.com
 JASCO Deutschland www.jasco.de | JASCO UK www.jasco.co.uk | JASCO France www.jascofrance.fr

JASCO Benelux www.jasco.nl | JASCO Spain www.jasco-spain.com

Algeria, Austria, Belgium, Cyprus, Denmark, Egypt, Finland, France, Germany, Greece, Hungary,
 Israel, Italy, Jordan, Kuwait, Luxembourg, Morocco, Netherlands, Norway, Poland, Portugal, Romania,
 Saudi Arabia, South Africa, Spain, Sweden Switzerland, Tunisia, Turkey,
 United Arab Emirates, United Kingdom, Yemen

JASCO CHINA (SHANGHAI) CO., LTD.

Room No.D, 10F, World Plaza, 855 Pudong South Road, Pudong New Area, Shanghai, China
 Tel: +86-21-6888-7871 Fax: +86-21-6888-7879 Web: www.jasco-global.com
 China



Products described herein are designed and manufactured by ISO-9001 and ISO-14001 certified JASCO Corporation